
- **Rotores disponibles**

Fixed angle rotor					
Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	24×1.5/2ml (aerosol tight)	0494021100	15200rpm/21697×g	24×1×0.2ml	1101010001
				24×1×0.5ml	1101010002
	12×5ml (aerosol tight)	0494030500	14500rpm/19980×g	12×1×1.8ml frozen storage tube	1101020002
				12×1×2ml frozen storaget ube	1101020004
				12×1×2ml	1101020001
	8×8×0.2ml PCR strip (aerosol tight)	0494010900	15200rpm/18856×g	/	/