- **AVAILABLE ROTORES**

Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	12×1.5/2ml	0494020100	16500rpm/18262×g	12×1×0.2ml	1101010001
				12×1×0.5ml	1101010002
	8×5ml	0494030100	16500rpm/18262×g	8×1×1.8ml frozen storage tube	1101020002
				8×1×2ml frozen storage tube	1101020004
				8×1×2ml	1101020001
	3×8×0.2ml PCR strip	0494010200	16500rpm/18262×g	/	/