




- **AVAILABLE ROTORS**

Fixed angle rotor					
Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	12×1.5/2ml	0494020300	20000rpm/26742×g	12×1×0.2ml	1101010001
				12×1×0.5ml	1101010002
	24×1.5/2ml aerosol tight	0495210100	16000rpm/24327×g	24×1×0.2ml	1101010001
				24×1×0.5ml	1101010002
	48×1.5/2ml aerosol tight	0494020900	13000rpm/18157×g	48×1×0.2ml	1101010001
				48×1×0.5ml	1101010002
	16×5ml round bottom (aerosol tight)	0494030600	15000rpm/21381×g	12×1×1.8ml frozen storage tube	1101020002
				12×1×2ml frozen storage tube	1101020004
				12×1×2ml	1101020001
	12×10ml round bottom (aerosol tight)	0494040600	14000rpm/17267×g	12×1×2ml	1101030001
				12×1×5ml	1101030002
	12×15ml conical bottom (aerosol tight)	0494050500	13000rpm/17212×g	12×1×5ml	1101040001
				12×1×10ml	1101040002
	6×50/15ml conical bottom (aerosol tight)	0494072200	12000rpm/17065×g	6×1×2ml	1101070003
				6×1×5ml	1101070004/conical tube 1101070005/round tube
				6×1×7ml	1101070012
				6×1×10ml	1101070006/glass tube 1101070007/conical tube 1101070008/round tube
				6×1×15ml	1101070001/conical tube
				6×1×25ml	1101070010
				6×1×50ml round tube	1101070002

• AVAILABLE ROTORS

Fixed angle rotor					
Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	8×50ml conical bottom (aerosol tight)	0494071400	11000rpm/13865×g	8×1×2ml	1101070003
				8×1×5ml	1101070004/conical tube 1101070005/round tube
				8×1×7ml	1101070012
				8×1×10ml	1101070006/glass tube 1101070007/conical tube 1101070008/round tube
				8×1×15ml	1101070001/conical tube
				8×1×25ml	1101070010
				8×1×50ml round tube	1101070002
	6×100ml round bottom (aerosol tight)	0494080600	10000rpm/11258×g	6×1×5ml	110108000
				6×1×10ml 6×1×15ml conical tube	7 110108000
				6×1×30ml	1101080002 1
					1101080003
				6×1×50ml	1101080005/conical tube 1101080004/round tube
	8×8×0.2ml PCR strip	0494010700	15000rpm/18363×g	/	/
	12x10ml blood tube	0494501900	5000rpm/3186×g	/	/

Swing out rotor					
Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	Capacity: 2x96 well plate	0495010600	4200rpm/1972×g	/	/
	Capacity: 24x1.5/2ml (aerosol tight)	0495210100	12000rpm/13845×g	/	/
	Capacity: 4x100ml	0495220100	4500rpm/3078×g	4×1×10ml	
				4x1x15ml conical bottom	
				4x1x50ml conical bottom	