




• AVAILABLE ROTORES

Fixed angle rotor					
Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	12×1.5/2ml	0494020200	19000rpm/ 24135×g	12×1×0.2ml	1101010001
				12×1×0.5ml	1101010002
	24×1.5/2ml	0494020600	16000rpm/24470×g	24×1×0.2ml	1101010001
				24×1×0.5ml	1101010002
	24×1.5/2ml aerosol tight	0494020500	17700rpm/29246×g	24×1×0.2ml	1101010001
				24×1×0.5ml	1101010002
	36×1.5/2ml	0494020800	14000rpm/18275×g	36×1×0.2ml	1101010001
				36×1×0.5ml	1101010002
	48×1.5/2ml	0494020900	13000rpm/17741×g	48×1×0.2ml	1101010001
				48×1×0.5ml	1101010002
	12×5ml round bottom	0494030200	17000rpm/24232×g	12×1×1.8ml frozen storage tube	1101020002
				12×1×2ml frozen storage tube	1101020004
				12×1×2ml	1101020001
	6×10ml round bottom	0494040100	16000rpm/17601×g	6×1×2ml	1101030001
				6×1×5ml	1101030002
	12×10ml round bottom	0494040200	14000rpm/17070×g	12×1×2ml	1101030001
				12×1×5ml	1101030002
	12×15ml conical bottom	0494050200	12000rpm/14730×g	12×1×5ml	1101040001
				12×1×10ml	1101040002

AVAILABLE ROTORES


Fixed angle rotor					
Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	6×50ml round bottom	0494070100	13000rpm/17628×g	6×1×2ml	1101060003
				6×1×5ml	1101060004
				6×1×10ml	1101060001
				6×1×15ml conical tube	1101060005
				6×1×30ml	1101060002
	6×50ml conical bottom	0494070200	13000rpm/18232×g	6×1×2ml	1101070003
				6×1×5ml	1101070004/conical tube 1101070005/round tube
				6×1×7ml	1101070012
				6×1×10ml	1101070006/glass tube 1101070007/conical tube 1101070008/round tube
				6×1×15ml conical tube	1101070001
				6×1×25ml	1101070010
				6×1×50ml round tube	1101070002
	8×50ml round bottom	0494070400	11000rpm/12513×g	8×1×2ml	1101060003
				8×1×5ml	1101060004
				8×1×10ml	1101060001
				8×1×15ml conical tube	1101060005
				8×1×30ml	1101060002
	8×50ml conical bottom	0494070500	11000rpm/13365×g	8×1×2ml	1101070003
				8×1×5ml	1101070004/conical tube 1101070005/round tube
				8×1×7ml	1101070012
				8×1×10ml	1101070006/glass tube 1101070007/conical tube 1101070008/round tube
				8×1×15ml conical tube	1101070001
				8×1×25ml	1101070010
				8×1×50ml round tube	1101070002

AVAILABLE ROTORES


Fixed angle rotor

Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	12×50ml conical bottom	0494070700	9000rpm/11057×g	12×1×2ml 12×1×5ml 12×1×7ml 12×1×10ml 12×1×15ml 12×1×25ml 12×1×50ml	1101070003 1101070004/conical tube 1101070005/round tube 1101070012 1101070006/glass tube 1101070007/conical tube 1101070008/round tube 1101070001/conical tube 1101070010 1101070002/round tube 1101080007
	4×100ml round bottom	0494080100	11000rpm/12743×g	4×1×5ml 4×1×10ml 4×1×15ml conical tube 4×1×30ml 4×1×50ml	1101080001 1101080002 1101080003 1101080005/conical tube 1101080004/round tube 1101080007 1101080001
	6×100ml round bottom	0494080200	10000rpm/11604×g	6×1×5ml 6×1×10ml 6×1×15ml conical tube 6×1×30ml 6×1×50ml	1101080002 1101080003 1101080005/conical tube 1101080004/round tube

Swing out rotor

Rotor	Capacity	Order No.	Max speed/RCF	Adaptor	
				Capacity	Order No.
	Capacity: 2×3×96 well (micro plate)	0495011000	4000rpm/2146×g	2×24×1.5/2ml	1102040001

Swing out rotor

Rotor body	Holder	Order No. (Rotorbody+holder)	Max speed/RCF	Adaptor	
				Capacity	Order No.
	 Capacity: 4×300ml Order No.:0495050220 Note: No suitable bottle,must choose adaptor	0495050200	4500rpm/3826×g	4×15×1.5/2ml 4x12x5ml 4×16×7ml blood tube 4x10x10ml 4x9x10ml blood tube 4×9×15ml conical tube 4×4×50ml conical tube 4×1×100ml 4×1×250ml	No. 0497012600 1102060011 0497011700 1102060012 1102060013 0497011200 0497011300 0497012500 1102060005
Order No.: 0495180510	 Capacity: 4x250ml conical bottle Order No.:1004010017	0495050200	4500rpm/4369×g	/	/

• AVAILABLE ROTORS

Swing out rotor					
Rotor body	Holder	Order No. (Rotorbody+holder)	Max speed/RCF	Adaptor	
				Capacity	Order No.
	 Capacity: 4×500ml Order No.: 1006010009 (without lid) Aerosol tight lid Order No.: 0495061110 Note: Aerosol tight lid cannot be used when use 500ml bottle	0495060400 (without lid) 0495061300 (with lid)	4200rpm/3551×g	4×18×2ml	0497012700
				4×18×2ml frozen storage tube	0497012900
				4×19×7ml blood tube	0497010200
				4×12×10ml	0497012800
				4×14×10ml blood tube	1102070007
				4×9×15ml conical tube	0497010300
				4×4×50ml conical tube	0497010400
				4×2×100ml	049701300
				4×1×100ml	0
				4×1×250ml	110207000
 Order No.: 0495180710	 Capacity: 4×500ml (square bucket) Order No.: 0495060510 Note: No suitable bottle,must choose adaptor	0495060500	4200rpm/3551×g	4×25×2ml	5
				4×25×2ml frozen storage tube	110207000 04970135006 049701320
				4×22×5ml	0497013300
				4×28×7ml blood tube	0497010100 0
				4×20×10ml blood tube	1102080003
				4×10×15ml conical tube	0497010500
				4×6×30ml	0497013400
				4×5×50ml conical tube	0497010600
				4×2×100ml	110208000
				4×1×250ml	6
	 Capacity: 4×2×96 well (micro plate) Order No.: 0495060610	0495060600	4200rpm/3551×g	4×24×1.5/2ml	110208000 7 1102040001
	 Capacity: 96×2~7ml (blood tube) Order No.: 0495180720	0495180700	4200rpm/3551×g	/	/
	 Capacity: 148×5ml (RIA tube) Order No.: 0496061010	0495061000	4200rpm/3551×g	/	/